



bio-chem – The optimal cleaners and degreasers for use in the BIO-CIRCLE Clean Box



- ➔ Degreasing with CB 100
- ➔ Removal of special greases and waxes with GS 200
- ➔ Removal of old paint and varnish with GS 200
- ➔ Neutral derusting with BIO-RUST
- ➔ Dewatering with OMNI
- ➔ Removal of white rust with Limescale Barracuda
- ➔ Flushing after rinsing
- ➔ Rinsing bath with corrosion protection for water-based systems

NEW BIO-CIRCLE[®] Clean Box Flow

Cleaning directly on the work piece.

The Clean Box Flow is a cleaning system that can be used flexibly. With its mobile trolley the Clean Box Flow can be placed directly next to the parts to be cleaned. The battery-powered pump and the washing brush facilitate easy cleaning in virtually any place. The cleaning medium is filtered during the cleaning of the work piece, which ensures longer service life.

The Clean Box Flow is sturdy and has been developed for practical use with the bio-chem cleaners. The mobile trolley features storing space for jerrycans or work pieces. The Clean Box Flow consists of the immersion box with a 12 V pump, storage battery and washing brush, the lid and the trolley.

- ➔ Manual cleaning with the washing brush
- ➔ Mobile trolley
- ➔ Ergonomic working height
- ➔ Storage battery and charger included
- ➔ Storing space for jerrycans or work pieces
- ➔ Battery-powered pump, easy cleaning
- ➔ Cleaning media: CB 100, CB 100 Alu, and CB 100 F



Clean Box Flow
– with integrated pump and storage battery

